**SERALENE®**

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**SUTURE MATERIALS
OF POLYVINYLIDENE FLUORIDE**

- SERALENE® blue

DESCRIPTION

SERALENE® is a synthetic monofilament nonabsorbable sterile surgical suture material. The basic material of the thread is polyvinylidene fluoride with the molecular formula $(\text{CH}_2-\text{CF}_2)_n$, which is synthesized from 1,1-difluoroethylene. SERALENE® blue is dyed with solvent blue, Colour Index No. 61568.

Our polyvinylidene fluoride suture materials are supplied as combinations with stainless steel needles. They are also supplied in a number of variants with specific accessories for intracutaneous sutures and tenorrhaphy. They satisfy the requirements of Directive 93/42/EEC on medical devices and the harmonized standard of the European Pharmacopoeia for "sutures, sterile nonabsorbable (fila non resorbilia sterilis)".

INDICATIONS

SERALENE® is suitable for adaptation of soft tissues and for ligatures where nonabsorbable suture material is indicated, including use in cardiovascular and neurological procedures as well as in microsurgery and ophthalmic procedures. Polyvinylidene fluoride threads can also be used as holding sutures and for marking purposes.

USE

In the choice and use of suture material, due account should be taken of the condition of the patient, the experience of the doctor, the surgical technique, and the size of the wound. Standard knot-tying techniques should be used to ensure knot-holding security. In order to prevent needle damage during suturing it is recommended that the needle be grasped at a point between the end of its first third (as seen from its end) and its midpoint.

ACTION

The function of a medical suture in wound care is to attach tissue to tissue or tissue to allogeneic material. Ligatures are a special type of suture used to close hollow organs.

After a slight initial tissue reaction, suture material of polyvinylidene fluoride is slowly and progressively encapsulated by connective tissue. It is nonabsorbable and physiologically inert and therefore retains its tensile strength indefinitely. Because of its monofilament thread structure, SERALENE® can also be used in tissue at risk of infection.

CONTRAINDICATIONS

None known.

WARNINGS / PRECAUTIONS / INTERACTIONS

The requirements of the suture material and the techniques needed vary with the indication. Therefore, the user should be familiar with the surgical techniques to be employed. In the choice of suture material, account should be taken of the in-vivo characteristics of the suture material concerned. Like all foreign bodies, SERALENE® can give rise to calculi when in prolonged contact with the biliary or urinary tract. In the case of sutures in infected tissue the use of suitable ancillary measures is recommended. Surgical instruments should be handled with care in order to avoid damage to threads. Bending of needles or grasping them outside of the recommended area can result in needle breakage. Because of the risk of infection from needle pricks, needles should be handled and disposed of with particular care.

ADVERSE REACTIONS

The following adverse reactions can occur in association with correct use of SERALENE®: transient local irritation, inflammatory reactions to the foreign body.

STERILITY

SERALENE® is supplied sterile (sterilization method: ethylene oxide). Do not resterilize, critical changes in suitability for use possible (more de-

tails available on request)! Do not use damaged packs! Discard unused opened packs!

STORAGE CONDITIONS

No special storage conditions are required. Do not use after the expiry date!

PRESENTATION

SERALENE® is supplied in various diameters and lengths in combination with various types of needle. The combinations are supplied in individual packs and multiple packs in a sterile package, and in the case of intracutaneous sutures and tenorrhaphy can be supplied with specific accessories. Pledgets to reinforce the body's own structures may be enclosed with SERALENE®. More precise details are given in the catalogue.

EXPLANATION OF SYMBOLS ON THE PACKAGE

 CE symbol and identification number of the notified body. The product satisfies the essential requirements of Directive 93/42/EEC on medical devices

 Surgical suture, nonabsorbable, monofilament, dyed

 Polyvinylidenfluoride

 Antireflective needle

 Detachable needle

 Pledget (length times width times thickness [mm])

Needle identification symbols and abbreviations are explained in the catalogue.

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SERAG-WIESSNER GmbH & Co. KG