

ÖHRNADELN *SURGICAL EYE NEEDLES*

Chirurgisches Nahtmaterial
Surgical Sutures



SERAG
WIESSNER

Die Bezeichnungen unserer Nadeln setzen sich aus einer Buchstaben-Zahlenkombination nach den Empfehlungen des Technischen Ausschusses der Vereinigung der Hersteller von chirurgischem Nahtmaterial e.V. zusammen. Der erste Buchstabe gibt die Nadelform an, der zweite Buchstabe kennzeichnet die Nadelart. Folgt ein dritter bzw. vierter Buchstabe, so bezieht sich dieser auf besondere Merkmale der Nadel. Die Zahl nach den Buchstaben gibt in mm die Nadellänge in gestreckter Form an.

The designations of our surgical needles consist of a letter-number-combination as per suggestion of the Technical Committee of the Association of Manufacturers of Surgical Sutures. The first letter indicates the needle shape, the second letter indicates the needle type. If a third or fourth letter follows, these refer to special characteristics of the needle. The number after the letters indicates the overall length of the needle in mm.

Nadelform · Needle shape

D =	$\frac{3}{8}$ kreisförmig	$\frac{3}{8}$ circle
G =	gerade Nadel	straight
H =	$\frac{1}{2}$ kreisförmig	$\frac{1}{2}$ circle

Nadelart · Needle type


R =	Rundkörpernadel	round-bodied
S =	schneidender Nadelkörper	reverse cutting

Besondere Merkmale · Special characteristics

F =	feine Nadel	fine needle
F	(nach den Zahlen)	(behind the numbers)
=	extra feine Nadel	extra fine needle
T =	Trokarnadel	trocarn needle
X =	extra stark	extra strong

Beispiele · Examples:

DS 18	D	$\frac{3}{8}$ kreisförmig	$\frac{3}{8}$ circle
	S	schneidender Nadelkörper	reverse cutting
	18	18 mm lang (gestreckte Länge)	18 mm long (straightened length)
HRX 22	H	$\frac{1}{2}$ kreisförmig	$\frac{1}{2}$ circle
	R	Rundkörpernadel	round-bodied
	X	extra starke Nadel	extra strong
	22	22 mm lang (gestreckte Länge)	22 mm long (straightened length)



Federöhrnadel surgical needle spring eye



Nadelbezeichnung Needle	Symbol des Nadelquerschnitts Symbol of section	Beschreibung Description
DS		$\frac{3}{8}$ kreisförmig, schneidender Nadelkörper $\frac{3}{8}$ circle, reverse cutting needle
DSF		$\frac{3}{8}$ kreisförmig, schneidender Nadelkörper (feine Nadel) $\frac{3}{8}$ circle, reverse cutting needle (fine needle)
GR		gerade, Rundkörpernadel straight, round-bodied needle
GS		gerade, schneidender Nadelkörper straight, cutting needle
HR		$\frac{1}{2}$ kreisförmig, Rundkörpernadel $\frac{1}{2}$ circle, round-bodied needle
HRF		$\frac{1}{2}$ kreisförmig, Rundkörpernadel (feine Nadel) $\frac{1}{2}$ circle, round-bodied needle (fine needle)
HR- ..-F		$\frac{1}{2}$ kreisförmig, Rundkörpernadel (extra feine Nadel) $\frac{1}{2}$ circle, round-bodied needle (extra fine needle)
HRT		$\frac{1}{2}$ kreisförmig, Rundkörpernadel mit Trokarspitze (Trokarnadel) $\frac{1}{2}$ circle, round-bodied needle with trocar point (trocar needle)
HRX		$\frac{1}{2}$ kreisförmig, Rundkörpernadel (extra starke Nadel) $\frac{1}{2}$ circle, round-bodied needle (extra strong needle)
HS		$\frac{1}{2}$ kreisförmig, schneidender Nadelkörper $\frac{1}{2}$ circle, reverse cutting needle
HSF		$\frac{1}{2}$ kreisförmig, schneidender Nadelkörper (feine Nadel) $\frac{1}{2}$ circle, reverse cutting needle (fine needle)
HSX		$\frac{1}{2}$ kreisförmig, schneidender Nadelkörper (extra starke Nadel) $\frac{1}{2}$ circle, reverse cutting needle (extra strong needle)

ÖHRNADELN / SPRING EYE

DS		Bestellnummer Code number steril sterile	Fadenzuordnung USP suture size recommendation USP	Fadenzuordnung EP suture size recommendation EP
	 DS-14	3114	5/0 - 3/0	1 - 2
	 DS-18	3113	5/0 - 3/0	1 - 2
	 DS-22	3112	5/0 - 3/0	1 - 2
	 DS-24	3111	5/0 - 3/0	1 - 2
	 DS-27	3110	5/0 - 3/0	1 - 2
	 DS-33	3109	5/0 - 3/0	1 - 2
	 DS-36	3108	2/0 - 0	3 - 3,5
	 DS-40	3107	2/0 - 0	3 - 3,5
	 DS-43	3106	0 - 2	3,5 - 5
	 DS-47	3105	0 - 2	3,5 - 5
	 DS-52	3104	0 - 2	3,5 - 5
	 DS-56	3103	0 - 2	3,5 - 5
	 DS-63	3102	0 - 2	3,5 - 5
	 DS-68	319S	1 - 3+4	4 - 6
	 DS-75	319R	1 - 3+4	4 - 6
	 DS-95	3101	1 - 3+4	4 - 6

Packungsinhalt: steril - 24 Folien zu jeweils 2 Stück (= 48 Stück)
 Packaging unit: sterile - 24 units with 2 needles each

DSF fein fine		Bestellnummer Code number steril sterile	Fadenzuordnung USP suture size recommendation USP	Fadenzuordnung EP suture size recommendation EP
	DSF-16 	3174	5/0 - 4/0	1 - 1,5
	DSF-21 	3173	5/0 - 3/0	1 - 2
	DSF-25 	3172	5/0 - 3/0	1 - 2
	DSF-29 	3171	5/0 - 3/0	1 - 2
	DSF-36 	3170	5/0 - 3/0	1 - 2





GS		Bestellnummer Code number steril sterile	Fadenzuordnung USP suture size recommendation USP	Fadenzuordnung EP suture size recommendation EP
	GS-42 	3189	0 - 2	3,5 - 5
	GS-45 	319E	0 - 2	3,5 - 5

GR		Bestellnummer Code number steril sterile	Fadenzuordnung USP suture size recommendation USP	Fadenzuordnung EP suture size recommendation EP
	GR-22 	3190	5/0 - 3/0	1 - 2
	GR-36 	3191	5/0 - 3/0	1 - 2
	GR-60 	319Q	2/0 - 0	3 - 3,5

ÖHRNADELN / SPRING EYE

Packungsinhalt: steril - 24 Folien zu jeweils 2 Stück (= 48 Stück)
Packaging unit: sterile - 24 units with 2 needles each

ÖHRNADELN / SPRING EYE

HR-..F extra fein extra fine		Bestellnummer Code number steril sterile	Fadenzuordnung USP suture size recommendation USP	Fadenzuordnung EP suture size recommendation EP
	HR-15-F	3187	5/0 - 3/0	1 - 2
	HR-17-F	3186	5/0 - 3/0	1 - 2
	HR-20-F	3185	5/0 - 3/0	1 - 2
	HR-23-F	3184	5/0 - 3/0	1 - 2
	HR-30-F	3183	5/0 - 3/0	1 - 2
	HR-33-F	3182	5/0 - 3/0	1 - 2
	HR-37-F	3181	5/0 - 3/0	1 - 2

HRF fein fine		Bestellnummer Code number steril sterile	Fadenzuordnung USP suture size recommendation USP	Fadenzuordnung EP suture size recommendation EP
	HRF-20	319F	5/0 - 3/0	1 - 2
	HRF-22	3124	5/0 - 3/0	1 - 2
	HRF-26	3123	2/0 - 0	3 - 3,5
	HRF-35	3122	2/0 - 0	3 - 3,5
	HRF-40	3121	2/0 - 0	3 - 3,5
	HRF-45	3120	2/0 - 1	3 - 4
	HRF-50	3119	2/0 - 1	3 - 4
	HRF-65	3118	2/0 - 1	3 - 4
	HRF-70	3117	2/0 - 1	3 - 4

HR



Bestellnummer
Code number
steril
sterile

Fadenzuordnung USP
suture size
recommendation USP

Fadenzuordnung EP
suture size
recommendation EP

HR-15

3146

5/0 - 3/0

1 - 2

HR-17

3145

5/0 - 3/0

1 - 2

HR-20

3144

5/0 - 3/0

1 - 2

HR-22

3143

5/0 - 3/0

1 - 2

HR-26

3142

2/0 - 0

3 - 3,5

HR-32

3141

2/0 - 0

3 - 3,5

HR-38

3140

2/0 - 0

3 - 3,5

HR-42

3139

0 - 2

3,5 - 5

HR-45

3138

0 - 2

3,5 - 5

HR-47

319N

0 - 2

3,5 - 5

HR-52

319M

0 - 2

3,5 - 5

HR-55

319L

0 - 2

3,5 - 5

HR-62

319K

1 - 3+4

4 - 6

HR-65

319J

1 - 3+4

4 - 6

HR-70

319H

1 - 3+4

4 - 6

HR-75

319G

1 - 3+4

4 - 6

HR-95

319Z

1 - 3+4

4 - 6

ÖHRNADELN / SPRING EYE

Packungsinhalt: steril - 24 Folien zu jeweils 2 Stück (= 48 Stück)
Packaging unit: sterile - 24 units with 2 needles each

ÖHRNADELN / SPRING EYE

HRT		Bestellnummer Code number steril sterile	Fadenzuordnung USP suture size recommendation USP	Fadenzuordnung EP suture size recommendation EP
	HRT-20	3169	2/0 - 0	3 - 3,5
	HRT-22	3168	2/0 - 0	3 - 3,5
	HRT-26	3167	2/0 - 0	3 - 3,5
	HRT-30	3166	0 - 2	3,5 - 5
	HRT-38	3165	0 - 2	3,5 - 5
	HRT-42	3164	1 - 3+4	4 - 6
	HRT-48	3163	1 - 3+4	4 - 6
	HRT-58	316Z	1 - 3+4	4 - 6
	HRT-64	316Y	1 - 3+4	4 - 6
	HRT-85	316W	1 - 3+4	4 - 6

Packungsinhalt: steril - 24 Folien zu jeweils 2 Stück (= 48 Stück)
Packaging unit: sterile - 24 units with 2 needles each



Bestellnummer
Code number
steril
sterile

Fadenzuordnung USP
suture size
recommendation USP

Fadenzuordnung EP
suture size
recommendation EP

HRX-20



319T

2/0 - 0

3 - 3,5

HRX-25



3188

0 - 2

3,5 - 5

HRX-27



3129

0 - 2

3,5 - 5

HRX-29



3128

0 - 2

3,5 - 5

HRX-34



319W

0 - 2

3,5 - 5

HRX-38



3127

0 - 2

3,5 - 5

HRX-40



3126

0 - 2

3,5 - 5

HRX-46



3125

1 - 3+4

4 - 6

HRX-48



319P

1 - 3+4

4 - 6

HRX-55



319U

1 - 3+4

4 - 6

HRX-68



319V

2 - 5

5 - 7



Bestellnummer
Code number
steril
sterile

Fadenzuordnung USP
suture size
recommendation USP

Fadenzuordnung EP
suture size
recommendation EP

HSF-12



3178

5/0 - 3/0

1 - 2

HSF-15



3177

5/0 - 3/0

1 - 2

HSF-17



3176

5/0 - 3/0

1 - 2

ÖHRNADELN / SPRING EYE

Packungsinhalt: steril - 24 Folien zu jeweils 2 Stück (= 48 Stück)
Packaging unit: sterile - 24 units with 2 needles each

ÖHRNADELN / SPRING EYE

HS		Bestellnummer Code number steril sterile	Fadenzuordnung USP suture size recommendation USP	Fadenzuordnung EP suture size recommendation EP
	HS-15	3162	5/0 - 3/0	1 - 2
	HS-17	3161	5/0 - 3/0	1 - 2
	HS-20	3160	5/0 - 3/0	1 - 2
	HS-23	3159	2/0 - 0	3 - 3,5
	HS-26	3158	2/0 - 0	3 - 3,5
	HS-30	3157	2/0 - 0	3 - 3,5
	HS-33	3156	2/0 - 0	3 - 3,5
	HS-36	3155	0 - 2	3,5 - 5
	HS-40	3154	0 - 2	3,5 - 5
	HS-48	3153	0 - 2	3,5 - 5
	HS-50	3152	0 - 2	3,5 - 5
	HS-55	3151	0 - 2	3,5 - 5
	HS-60	3150	1 - 3+4	4 - 6
	HS-66	3149	1 - 3+4	4 - 6
	HS-70	3148	1 - 3+4	4 - 6
	HS-75	3147	1 - 3+4	4 - 6

Packungsinhalt: steril - 24 Folien zu jeweils 2 Stück (= 48 Stück)
 Packaging unit: sterile - 24 units with 2 needles each



Bestellnummer Code number steril sterile	Fadenzuordnung USP suture size recommendation USP	Fadenzuordnung EP suture size recommendation EP
---	---	---

HSX-20



319D

2/0 - 0

3 - 3,5

HSX-25



3137

0 - 2

3,5 - 5

HSX-27



3136

0 - 2

3,5 - 5

HSX-29



3135

0 - 2

3,5 - 5

HSX-34



319X

0 - 2

3,5 - 5

HSX-38



3134

0 - 2

3,5 - 5

HSX-40



3133

1 - 3+4

4 - 6

HSX-46



3132

1 - 3+4

4 - 6

HSX-48



3131

1 - 3+4

4 - 6

HSX-55




3130

2 - 5

5 - 7

HSX-60




319C

2 - 5

5 - 7

HSX-70




319B

2 - 5

5 - 7

HSX-75



319A

3+4 - 8

6 - 10

ÖHRNADELN / SPRING EYE

Packungsinhalt: steril - 24 Folien zu jeweils 2 Stück (= 48 Stück)
Packaging unit: sterile - 24 units with 2 needles each



SERAG-WIESSNER GmbH & Co. KG
Zum Kugelfang 8 - 12
95119 Naila/Germany

 + 49 9282 937-0

 + 49 9282 937-9369

Export Department:

 + 49 9282 937-230

 + 49 9282 937-9785

 info@serag-wiessner.de

 www.serag-wiessner.de