

Chirurgisches Nahtmaterial *Surgical Sutures*

SERAFAST®

PGACL

POLYGLYKOLSÄURE-CAPROLACTON
POLYGLYCOLIC ACID-CAPROLACTONE



violett / monofil / resorbierbar
violet / monofilament / absorbable















ungefärbt / monofil / resorbierbar
undyed / monofilament / absorbable

Resorptionsprofil: 50 % Reißkraft nach ca. 8 - 13 Tagen
Abbauvorgang abgeschlossen nach ca. 90 - 120 Tagen

*Absorption profile: 50 % tensile strength after 8 - 13 days
Absorption process completed after
approx. 90 - 120 days*












**SERAG
WIESSNER**

Nadelbezeichnung Needle		Stärke USP Size USP	Stärke metric Size metric	Fadenlänge/Meter Length of suture/metre	Abgabemenge/Stück Units per box	Bestell-Nummer ungefärbt Code number undyed	Bestell-Nummer violett Code number violet
DS-12		5/0	1	1 x 0,45	24	40101312	
DS-15		5/0	1	1 x 0,70	24		JO103213
		4/0	1,5	1 x 0,70	24	40153213	
DS-18		5/0	1	1 x 0,45	24	40101314	
		5/0	1	1 x 0,70	24	40103214	JO103214
		4/0	1,5	1 x 0,45	24	40151314	
		4/0	1,5	1 x 0,70	24	40153214	JO153214
		3/0	2	1 x 0,45	24	40201314	
		3/0	2	1 x 0,70	24	40203214	JO203214
DS-25		4/0	1,5	1 x 0,70	24	40153215	JO153215
		3/0	2	1 x 0,70	24	40203215	JO203215
		2/0	3	1 x 0,70	24	40303215	JO303215
DS-30		2/0	3	1 x 0,70	24	40303216	
DSS-10		6/0	0,7	1 x 0,45	24	4007131S	
DSS-10 +DSS-10		6/0	0,7	1 x 0,70	24	4P07321S	
DSS-13		6/0	0,7	1 x 0,45	24	4007131L	
		5/0	1	1 x 0,45	24	4010131L	JO10131L
		5/0	1	1 x 0,70	24	4010321L	
		4/0	1,5	1 x 0,45	24	4015131L	
DSS-18		6/0	0,7	1 x 0,45	24	4007131M	
		4/0	1,5	1 x 0,45	24	4015131M	
		4/0	1,5	1 x 0,70	24	4015321M	
		3/0	2	1 x 0,45	24	4020131M	
		3/0	2	1 x 0,70	24	4020321M	
DSS-24		4/0	1,5	1 x 0,45	24	4015131N	
		3/0	2	1 x 0,45	24	4020131N	
		3/0	2	1 x 0,70	24	4020321N	
FRX-27 +FRX-27		3/0	2	1 x 0,70	24	4P203226	

SERAFAST®

NADEL-FADEN-KOMBINATIONEN
NEEDLE SUTURE COMBINATIONS


Nadelbezeichnung Needle		Stärke USP Size USP	Stärke metric Size metric	Fadenlänge/Meter Length of suture /metre	Abgabemenge/Stück Units per box	Bestell-Nummer ungefärbt Code number undyed	Bestell-Nummer violett Code number violet
GR-20		4/0	1,5	1 x 0,70	24	40153229	
GS-12		4/0	1,5	1 x 0,45	24	40151333	
GS-60		3/0	2	1 x 0,70	12	40203236	
HR-10		6/0	0,7	1 x 0,70	24	40073237	
HR-17		5/0	1	1 x 0,45	24	40101376	
		5/0	1	1 x 0,70	24	40103276	JO103276
		4/0	1,5	1 x 0,70	24	40153276	JO153276
		3/0	2	1 x 0,70	24	40203276	JO203276
HR-22		4/0	1,5	1 x 0,70	24	4015322C	JO15322C
		3/0	2	1 x 0,70	24	4020322C	JO20322C
HR-22 +HR-22		4/0	1,5	1 x 0,70	24		JP15322C
HR-26		4/0	1,5	1 x 0,70	24	40153241	JO153241
		3/0	2	1 x 0,70	24	40203241	JO203241
		2/0	3	1 x 0,70	24	40303241	JO303241
HR-36		2/0	3	1 x 0,70	24	40303243	


SERAFAST®

NADEL-FADEN-KOMBINATIONEN
NEEDLE SUTURE COMBINATIONS





SERAG-WIESSNER GmbH & Co. KG
Zum Kugelfang 8 - 12
95119 Naila/Germany

 + 49 9282 937-0

 + 49 9282 937-9369

Export Department:

 + 49 9282 937-230

 + 49 9282 937-9785

 info@serag-wiessner.de

 www.serag-wiessner.de